

Dark Dry Eschar below the knee

Except heels



Non-viable (dead) wound tissue that is characterized by a dry, thick leathery, black crust.

Interventions:

- Check **ABIs** first:
- If ABIs below 0.59 paint with Betadine BID
- If ABIs above 0.6: Autolytic or surgical removal of non-viable tissue

*Surgical only by healthcare professional

	Medicaid	99¢
Autolytic	Medihoney paste daily or Triad Q3 days cover with occlusive dressing	Zinc Oxide cream mixed with Petroleum Jelly daily cover with occlusive dressing

Wash wound with soap and warm water (preferably in shower) and pat dry between dressing changes

Product information

Medical Formulary

99 ¢

Medihoney or Triad

- Controls odor
- Debrides dry necrotic wounds with osmotic pull
- Dries out moist ischemic wounds



OR



Zinc Oxide mixed with Petroleum Jelly

- Debrides dry necrotic wounds by softening dry tissue



AND



- How to use: apply directly to the eschar and cover with occlusive dressing

- How to use: apply directly to the eschar and cover with occlusive dressing

- How to order:
shieldhealthcare.com
ccsmed.com
byramhealthcare.com